



COVID-19 PCR Test

Why Accula™?

The Accula™ SARS-CoV-2 Test is a rapid polymerase chain reaction (PCR) and lateral flow test for the qualitative, visual detection of the coronavirus SARS-CoV-2 viral RNA. This test uses comfortable lower nasal swab collection. The Accula™ SARS-CoV-2 Test is manufactured by Mesa Biotech and locally distributed by NHD, LLC.

FDA EUA Authorized

Rapid Results

Compatible with Accula or Solaris Dock

Features

- Identify presence of COVID-19 RNA
- Reliable results in 30 minutes
- Anterior nasal specimen collection
- Storage conditions: 1~30°C
- Analyzer docks sold separately
- CPT Code: 87635
- LOINC Code: 95409-9



Kit Contents

- 25 test devices
- 25 buffer vials
- 25 sterile nasal swabs
- 2 pipettes
- High positive control
- Low positive control
- Negative control
- Package insert
- Quick reference instructions

Performance Measures

- 95.8% Sensitivity
- 100% Specificity

Available Now

We offer reliable production and delivery schedules. Our team is committed to providing the highest quality service at the best possible prices.



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc



COVID-19 PCR Test

Product Specifications and Dimensions

Kit: 25 tests / 3.4 lbs. / 11.75" x 11.25" x 4.125"

Case: 150 tests / 6 kits / 23.2 lbs. / 24" x 12.25" x 13"

Pallet (SM): 2600 tests / 104 kits / 577 lbs. / 48"x 40"x 58"

Pallet (LG.): 4500 tests / 180 kits / 721 lbs. / 50"x 40"x 71"

Attachments and Documentation (click below)

FDA EUA Letter

FDA Instructions for Use

Healthcare Provider Fact Sheet

Patient Fact Sheet

Ordering Process

CALL (512) 651-1977 or EMAIL sales@nhd.llc for market pricing.

Request a purchase order online at nhd.llc or from your sales representative.

PLEASE NOTE: Due to the ongoing and rapidly changing pandemic conditions, pricing is variable and delivery times cannot be guaranteed by FedEx or UPS. We appreciate your business and strive to provide the most reliable product inventory and order processing.



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc