



CorDx COVID-19 Ag Test | At Home OTC

The CorDx COVID-19 Ag Test is a lateral flow immunoassay device intended for the qualitative detection of nucleocapsid protein antigen from the SARS-CoV-2 virus, delivering results in just 10 minutes.

Reliable test results delivered in just 10 minutes

Industry-Leading Performance High Affordability Rapid Results High Performance: 89.8% sensitivity 99.8% specificity 100% agreement in PCR testing during clinical trials* Approved 18-month shelf-life

In-house development allows for maximum cost savings and affordability

This test is authorized for non-prescription home use with self-collected anterior nasal (nares) swab samples from individuals aged 14 years or older, or adult collected anterior nasal (nares) swab samples from individuals aged two years or older. This test is authorized for individuals with symptoms of COVID-19 within the first 7 days of symptom onset when tested at least twice over three days with at least 48 hours between tests, and for individuals without symptoms or other epidemiological reasons to suspect COVID-19, when tested at least three times over five days with at least 48 hours between tests.



Call 512-651-1977 | Email sales@nhd.llc | Visit www.NHDmedical.com
Distributed by NHD, LLC. 5321 Industrial Oaks Blvd., Austin, TX. 78735-8823



CorDx COVID-19 Ag Test | At Home OTC



Rapid Results

Reliable test results delivered in just **10 MINUTES**



Industry-Leading Performance

High Performance:
89.8% SENSITIVITY
99.8% SPECIFICITY
100% agreement in PCR testing during clinical trials
Approved
18-MONTH SHELF-LIFE



High Affordability

In-house development allows for maximum **COST SAVINGS AND AFFORDABILITY**

Product Specifications and Dimensions

	Qty	Weight (lbs)	Dimensions (in)
Kit	2	0.09	6.5 x 3 x .71
Case	512	24.03	23.2 x 13.6 x 12.6
Pallet	15,360	790	48 x 40 x 80

Attachments and Documentation

(click below)

Package Insert

Fact Sheet for Healthcare Professionals

Instructions For Use



Call 512-651-1977 | Email sales@nhd.llc | Visit www.NHDmedical.com
Distributed by NHD, LLC. 5321 Industrial Oaks Blvd., Austin, TX. 78735-8823