



INDICAID™ COVID-19 Rapid Antigen Test

Why INDICAID™?

INDICAID™ COVID-19 Antigen Tests are a rapid testing device that detects symptomatic and asymptomatic COVID-19 in just 20 minutes. With anterior nasal swabs, specimen collection is quick and comfortable in point-of-care settings and does not require a machine reader to interpret results. INDICAID™ is manufactured by Phase Scientific in China and locally distributed by NHD, LLC.

Features

- Identify current COVID-19 infection status
- FDA Emergency Use Authorized
- Reliable results in 20 minutes
- Simple nasal specimen collection
- Qualitative detection of SARS-CoV-2 antigens
- All-in-one package, including collection swabs
- CPT Code: 86328-QW



Kit Contents

- 25 test devices
- 25 buffer solution vials
- 25 sterile nasal swabs
- Instructions for use
- Quick reference guide

Performance Measures

- 84.4% Sensitivity
- 97% Specificity



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc



INDICAID™ COVID-19 Rapid Antigen Test

Product Specifications and Dimensions

Kit: 25 tests / <1 lb. / 9" x 6.55" x 3.25"

Case: 450 tests / 18 kits / 17 lbs. / 20.5" x 18.8" x 10"

Pallet: 9,000 tests / 360 kits / 380 lbs. / 41" x 39.5" x 57.5"

Attachments and Documentation (click below)

FDA EUA Letter

Instructions for Use + Quick Reference Guide

Technical Brief

Patient Fact Sheet

Healthcare Provider Fact Sheet



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc