



QuickVue® SARS Antigen TEST

Why QuickVue®?

QuickVue® SARS Antigen tests do not require a machine reader and take 10 minutes to get results. This test is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation. It test uses a gentle self-collected anterior nasal (nares) swab sample to determine a positive or negative COVID-19 result. The swab is swirled in a tube of reagent solution, then removed, before a test strip is inserted.

Available Now

QuickVue® SARS Antigen Tests are manufactured by Quidel® and locally distributed by NHD, LLC. We offer reliable production and delivery schedules. Our team is committed to providing the highest quality service at the best possible prices.



Features


- Identify your current COVID-19 infection status
- Reliable results in 10 minutes
- Direct anterior nare specimen collection
- Fast and easy to use in point of care settings
- Qualitative detection of SARS-CoV-2 antigens
- Storage conditions: 15-30°C
- All-in-one package, including collection swabs

Kit Contents

- 25 test devices
- 25 sterile anterior nasal swabs
- Positive and negative control swabs
- Package insert with simple instructions

Performance Measures

- 96.6% Positive Percent Accuracy
- 99.3% Negative Percent Accuracy



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc



QuickVue[®]

SARS Antigen TEST

Attachments and Documentation
(click below)



QuickVue[®]
SARS Antigen TEST
FDA Instructions for Use
HCP Fact Sheet
Patient Fact Sheet

Ordering Process

CALL (512) 651-1977 or EMAIL sales@nhd.llc for market pricing.

Request a purchase order online at nhd.llc or from your sales representative.

PLEASE NOTE: Due to the ongoing and rapidly changing pandemic conditions, pricing is variable and delivery times cannot be guaranteed by FedEx or UPS. We appreciate your business and strive to provide the most reliable product inventory and order processing.



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc