



Sienna COVID-19 Antigen Rapid Test


Why choose Sienna?

The Sienna COVID-19 Antigen Rapid Test is a nasopharyngeal swab test kit intended for use in point of care settings. Easy-to-read results appear in just five minutes on the testing strip, indicating current infection status with 98.9% specificity. Sienna-Clarity COVID-19 Antigen Rapid Tests are manufactured in Finland and locally distributed in the U.S. by NHD, LLC.

FDA EUA Authorized
Easy to Use
CLIA Waived

Features

- Detect presence of SARS-CoV-2 nucleocapsid protein
- Reliable results in 5 minutes
- Nasopharyngeal specimen collection
- Fast and easy to use in point-of-care settings
- Storage conditions: 2-8°C
- All-in-one package, including collection swabs
- CPT Code: 86328-QW



Kit Contents

- 25 test devices
- 25 sterile nasopharyngeal swabs
- 25 buffer tubes
- 1 positive and 1 negative control swab
- 1 Work Station
- Package insert
- Quick start guide

Performance Measures

- 87.5% Sensitivity
- >98.9% Specificity

Available Now

We offer reliable production and delivery schedules. Our team is committed to providing the highest quality service at the best possible prices.



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc



Sienna COVID-19 Antigen Rapid Test

Product Specifications and Dimensions

Kit: 25 tests / 8.5" x 5.3" x 3.3" / .86 lbs.

Case: 1,225 tests / 49 kits / 23.6" x 15.7" x 21.7" / 42 lbs.

SM Pallet: 9,800 tests / 8 cases / 392 kits / 47.2" x 31.5" x 47.2" / 336 lbs

LG Pallet: 14,700 tests / 12 cases / 588 kits / 47.2" x 31.5" x 71" / 503 lbs.

Attachments and Documentation (click below)

FDA EUA Letter

FDA Letter of Condition

Instructions for Use

Healthcare Provider Fact Sheet

Patient Fact Sheet

Ordering Process

CALL (512) 651-1977 or EMAIL sales@nhd.llc for market pricing.

Request a purchase order form online at nhd.llc or from your sales representative.

PLEASE NOTE: Due to the ongoing and rapidly changing pandemic conditions, pricing is variable and delivery times cannot be guaranteed by FedEx or UPS. We appreciate your business and strive to provide the most reliable product inventory and order processing.



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc