



Status Flu A&B Test Kit

Why Status Flu A&B?

The Status Flu A&B Test is an in vitro rapid qualitative test that detects influenza type A and type B antigens directly from nasal swab and nasopharyngeal swab. The Status Flu A&B Test is manufactured by Lifesign in the U.S.A. and locally distributed by NHD, LLC.



Features

- CLIA Waived for swab specimens
- FDA Cleared
- Differentiates between Influenza A and B
- Nasal or nasopharyngeal specimen
- Flu A CPT Code: 87804QW
- Flu B CPT Code: 87804QW-59

Kit Contents

- 25 Test Devices
- 25 Flocked Nasal Swabs
- 25 Pre-measured Developer Solution Capsules
- Positive and Negative Control Swabs

Attachments and Documentation (click below)

[Procedure Card](#)

[Package Insert](#)

[Package Insert Addendum](#)

[Safety Data Sheet](#)

Performance Measures

- Flu A Sensitivity: 89.2%
- Flu A Specificity: 99.4%
- Flu B Sensitivity: 86.4%
- Flu B Specificity: 99.0%

Available Now

We offer reliable production and delivery schedules. Our team is committed to providing the highest quality service at the best possible prices.



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc