



Status Mono Whole Blood Rapid Test

Why Status Mono?

The Status Mono Test is an immunoassay for the qualitative detection of infectious mononucleosis IgM heterophile antibodies in serum, plasma or whole blood. The Status Mono Test is manufactured by Lifesign in the U.S.A. and locally distributed by NHD, LLC.



Attachments and Documentation (click below)

[Whole Blood Package Insert](#)

[25pk Package Insert](#)

[30pk Package Insert](#)

[Safety Data Sheet](#)

Features

- CLIA Waived for whole blood
- FDA Cleared with no age restrictions
- Results in 8 minutes
- Single reagent procedure
- Whole blood CPT Code: 86308QW

Kit Contents

- 10, 25, or 30 Tests
- Controls included for 25-30 tests
- Plastic collection pipettes
- Single Reagent
- Package Insert

Performance Measures

- Sensitivity: 99.9%
- Specificity: 98.8%

Available Now

We offer reliable production and delivery schedules. Our team is committed to providing the highest quality service at the best possible prices.



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc