




Status™ COVID-19/Flu A&B Antigen Test

Why Status™?

The Status™ COVID-19/Flu A&B Test is a rapid immunoassay for the simultaneous direct detection and differential diagnosis of SARS-CoV-2, Influenza Type A and Type B Antigen from nasopharyngeal swab specimens. This test provides 3 results from just one test and helps differentiate diagnoses from similar symptom presentation. Status™ COVID-19/Flu A&B Tests are manufactured by lifeSign and locally distributed by NHD, LLC.



Features

- CLIA Waived
- Determine COVID-19 and flu A&B infection status
- Visual results in 15 minutes
- Nasopharyngeal specimen collection
- Storage conditions: 2~30°C

FDA EUA Authorized | Rapid Results | Easy to Use

Kit Contents

- 25 test devices
- 25 sterile nasopharyngeal swabs
- Package insert
- Quick reference instructions

Performance Measures

- COVID-19 Sensitivity: 93.9% / Specificity: 100%
- Flu A Sensitivity: 91.4% / Specificity: 95.7%
- Flu B Sensitivity: 87.6% / Specificity: 95.9%

Available Now

We offer reliable production and delivery schedules.
Our team is committed to providing the highest quality service at the best possible prices.



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc



StatusTM COVID-19/Flu A&B Antigen Test

Product Specifications and Dimensions

Kit: 25 tests / 1.1 lbs. / 11" x 6.5" x 3.5"

Case: 500 tests / 15 kits / 23 lbs. / 21" x 17" x 17"

Pallet: 9,000 tests / 18 cases / 270 kits / 415 lbs. / 48" x 40" x 52"

COVID-19 CPT: 87811QW

FLU A CPT: 87804QW

FLU B CPT: 87804QW-59

Attachments and Documentation (click below)

StatusTM

FDA EUA Letter

FDA Instructions for Use

Health Care Provider Fact Sheets

Patient Fact Sheets



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc