



Why Syphilis Health Check™?

The Syphilis Health Check™ Rapid Antibody Test is a rapid qualitative screening for human TP antibodies in Whole Blood, Serum or Plasma. The Syphilis Health Check™ Rapid Antibody Test is locally distributed by NHD, LLC.



Features

- First and only CLIA Waived rapid treponemal syphilis test for whole blood
- Visual results in 10 minutes
- Whole blood specimen collection
- Reimbursable CPT Code: 86780QW

Performance Measures

- >97% Accuracy

Kit Contents

- 20 Test Devices
- 20 Fixed Volume Pipettes
- 5 mL diluent - saline buffer, detergent and sodium azide (NaN₃ < 0.1%)
- Instruction Leaflet
- Quick Reference Guide

Attachments and Documentation

(click below)

[Package Insert](#) | [Quick Reference Instructions](#) | [MSDS](#)



CALL : 512-651-1977 | EMAIL: sales@nhd.llc

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | www.nhd.llc